

Methods

The ASCproject

The Angiosarcoma Project (ASCproject) is a patient-driven genomics research initiative that enables angiosarcoma patients to enroll through an online portal. Patients can elect to provide consent for acquisition of a saliva sample, a blood sample, medical records, and archived FFPE tissue samples. The data generated as a result are used to create a clinically-annotated genomic dataset that can be shared across multiple platforms. The following describes the methods used to generate the current dataset from The Angiosarcoma Project (48 samples, 36 patients) in September of 2018. Please refer to ascproject.org/data-release for links to data on public platforms.

Methods

Design

The Angiosarcoma Project website (ASCproject.org) enables patients to enroll in the study remotely. The website and all associated messaging were developed through iterative feedback from angiosarcoma patients. Patients registered with first name, last name, email address, and confirmation of a angiosarcoma diagnosis. Registrants were asked to complete a 17-question survey, with all questions optional, about their experiences with angiosarcoma (see link on ascproject.org/data-release). Registrants acknowledged that responses would be stored in a secure database. Included in this acknowledgement is an understanding that patients may be re-contacted and that they can withdraw their information as indicated in the pre-submission statement:

"I understand that the information I entered here will be stored in a secure database and may be used to match me to one or more research studies conducted by the Angiosarcoma Project. If the information that I entered matches a study being conducted by the Angiosarcoma Project, either now or in the future, I agree to be contacted about possibly participating. I understand that if I would like my information deleted from the database, now or in the future, I can email info@ascproject.org and my information will be removed from the database."

All patients who at any time submitted a written request to withdraw from the study were immediately exited by study staff.

A link to an electronic informed consent document for formal enrollment in the study (see link on ascproject.org/data-release) was sent to registrants who completed the survey. Registrants could direct any questions to study staff throughout the enrollment process. Email reminders were sent weekly for three weeks, and again at six weeks, to registrants who had not completed the consent process. Registrants who provided informed consent were then asked to complete a medical release form in order to provide their contact information, as well as a list of physicians and hospitals that provided their care for angiosarcoma. Email reminders were sent weekly for three weeks, and again at six weeks, to registrants who had not completed the medical release form. Upon completion, signed copies of consent and medical release forms were sent electronically to the registrants. All data collected online were stored in a secure database.

Saliva kits were mailed to registrants who provided informed consent and who lived in The United States or

Canada. Saliva kits were labeled with a unique two-dimensional barcode and a pre-paid Business Reply Label addressed to the Broad Institute Genomics Platform prior to shipment. Each barcode identifier was assigned to a participant prior to shipment. Participants provided a saliva sample by following the included instructions (see link on ascproject.org/data-release) and returned the kits free of charge. Saliva kits received at the Broad Institute were logged by their unique barcodes, and stored at room temperature until they were advanced to Whole Exome Sequencing (WES).

Blood kits were mailed to registrants who provided informed consent, opted in to the blood biopsy component of the study, and lived in The United States or Canada. Blood kits were labeled with a unique two-dimensional barcode and included a pre-paid FedEx ClinPak envelope addressed to the Broad Institute Genomics Platform prior to shipment. Each barcode identifier was assigned to a participant prior to shipment. Participants provided a blood sample by bringing the kit and instructions (see link on ascproject.org/data-release) with them to their next regularly scheduled clinical appointment and requesting a courtesy draw. If a courtesy draw was not possible, patients were given the option to go to Quest Diagnostics with a voucher for their complementary blood draw. Blood kits were returned free of charge. Blood samples received at the Broad Institute were logged by their unique barcodes and fractionated. The resultant plasma and buffy coats were stored at -80C and plasma samples were advanced to Ultra-low Pass Whole Genome Sequencing (ULP-WGS) to ascertain the presence of circulating tumor DNA (ctDNA). Select plasma samples were advanced to WES. Buffy coats were advanced to WES where no saliva sample was available to generate matched germline data.

Study staff called the hospitals and physicians' offices listed in each participant's medical release form to confirm the fax number for the medical records department. A detailed request for medical records including clinic notes from treating providers, angiosarcoma treatment data (including radiation and chemotherapy), pathology reports, operative reports, referrals, MD to MD exchange, and genetic testing reports from the date of diagnosis through the date the request was faxed to each facility. Medical records were received by fax, mail, or secure electronic message. All medical records were scanned and uploaded to a secure drive to facilitate abstraction. The date of each procedure (e.g. biopsy, resection, mastectomy, etc.), type of procedure, histology, and facility which performed the procedure were abstracted from all pathology reports in order to prioritize tissue samples for request for patients who opted in to the tissue component of the study. Angiosarcoma samples were flagged for request. Study staff were informed not to request the most recent biopsy in order to avoid exhausting samples that may be needed for future clinical care.

Study staff called the pathology departments associated with each tissue sample to confirm the fax number for tissue requests. A form requesting a minimum of 8 and maximum of 23 5-micron unstained slides, or one tissue block, and one Hematoxylin and Eosin stain (H&E) slide was faxed to each pathology department. Requests explicitly stated that no sample should be exhausted in order to fulfill the request.

Tissue samples were received at the Broad Institute by mail. H&E slides and three additional unstained slides were sent for expert pathology review to confirm the presence of angiosarcoma in each sample. Tissue samples received as blocks were labeled with unique numerical identifiers and sent to the Dana-Farber/Harvard Cancer Center Specialized Histopathology Services (SHS) Core to be cut into three 30-micron scrolls per block. These scrolls were then labeled with unique barcode identifiers and submitted to the Broad Institute Genomics Platform for WES and Transcriptome Capture (RNA Seq). Tissue samples received as unstained slides were logged, labeled with unique barcode identifiers, and submitted to the

Medical Record Abstraction

A data dictionary (see Appendix 1, pg. 8) comprising 7 fields based on expert pathologist review and 6 fields from medical records was developed. The date of primary diagnosis with angiosarcoma was defined as the date of first definitive confirmation of disease by biopsy. Dates were abstracted to the greatest level of detail available in the record. Dates reported in the medical record only as a month and year were abstracted as the first of the month while dates reported only as a year were abstracted as the first of January of that year. For all other fields, only data explicitly reported in the medical record were abstracted; no data were inferred. In order to protect patient confidentiality, all dates were reported in shared data sets as elapsed time relative to the date of primary diagnosis and ages were grouped into 5-year increments.

Patient-Reported Data

Patient-reported data (PRD) comes from information provided by patients in the 17-question intake survey (see link on ascproject.org/data-release). When applicable, data were cleaned to standardize format as well as protect patient confidentiality. In order to protect confidentiality, race categories that were reported fewer than five times in the dataset were reclassified as “other”.

Biological Sample Processing

Methods for DNA Isolation from Saliva

DNA was extracted via the Chemagic MSM I with the Chemagic DNA Blood Kit-96 from Perkin Elmer. This kit combines a chemical and mechanical lysis with magnetic bead-based purification.

Saliva samples were incubated at 50°C for 2 hours. The saliva was then transferred to a deep well plate placed on the Chemagic MSM I. The following steps were automated on the MSM I.

M-PVA Magnetic Beads were added to the saliva. Lysis buffer was added to the solution and mixed. The bead-bound DNA was then removed from solution via a 96-rod magnetic head and washed in three Ethanol-based wash buffers. The beads were then washed in a final water wash buffer. Finally, the beads were dipped in elution buffer to resuspend the DNA sample in solution. The beads were then removed from solution, leaving purified DNA eluate.

DNA samples were quantified using a fluorescence-based PicoGreen assay.

Methods for DNA Isolation from Whole Blood

DNA was extracted via the Chemagic MSM I with the Chemagic DNA Blood Kit-96 from Perkin Elmer. This kit combines a chemical and mechanical lysis with magnetic bead-based purification.

Whole Blood samples were incubated at 37°C for 5-10 minutes to thaw. The blood was then transferred to a deep well plate with protease and placed on the Chemagic MSM I. The following steps were automated on the MSM I.

M-PVA Magnetic Beads were added to the blood, protease solution. Lysis buffer was added to the solution and vortexed to mix. The bead-bound DNA was then removed from solution via a 96-rod magnetic head and washed in three Ethanol-based wash buffers to eliminate cell debris and protein residue. The beads were

then washed in a final water wash buffer. Finally, the beads were dipped in elution buffer to re-suspend the DNA. The beads were then removed from solution, leaving purified DNA eluate.

DNA samples were quantified using a fluorescence-based PicoGreen assay.

cfDNA Extraction from Whole Blood

Whole blood was collected in EDTA, CellSave, or Streck tubes and processed for plasma fractionation within 3 hours of blood draw. Blood tubes were centrifuged at 1900 g for 10 minutes and plasma was transferred to second tube before further centrifugation at 15000 g for 10 minutes. Supernatant plasma was stored at -80C until cfDNA extraction. cfDNA was extracted using the QIAasympohy DSP Circulating DNA Kit according to the manufacturer's instructions, with 6.3 mL of plasma as input and with a 60 uL DNA elution (Qiagen, 2017).

Library Construction

Library construction was performed using the KAPA Hyper Prep Kit according to the manufacturer's instructions (Kapa, 2016), with the following modifications: initial DNA input was normalized to 20 ng in 50 uL of TE buffer (10mM Tris HCl 1mM EDTA, pH 8.0) according to picogreen quantification. For adapter ligation, Illumina paired end adapters were replaced with palindromic forked adapters, purchased from Integrated DNA Technologies, with unique dual-indexed molecular barcode sequences to facilitate downstream pooling. Adapters were diluted to 3.75 uM before addition to the samples. During PCR, 10 cycles were used. During the post-enrichment SPRI cleanup, elution volume was reduced to 30uL to maximize library concentration.

Post Library Construction Quantification and Normalization

Library quantification was performed using the Invitrogen Quant-It broad range dsDNA quantification assay kit (Thermo Scientific Catalog: Q33130) with a 1:200 PicoGreen dilution. Following quantification, each library is normalized to a concentration of 25 ng/uL, using a 1X Low TE pH 7.0 solution.

Library Pool Creation for Ultra-low Pass Sequencing

In preparation for the sequencing of the ultra-low pass libraries (ULP), approximately, 2 uL of the normalized library is transferred into a new receptacle and further normalized to a concentration of 3ng/uL, again using a 1X Low TE pH 7.0 solution. Following normalization, up to 96 individual ultra-low pass WGS pool is created via equivolume pooling.

In-solution hybrid selection for exome or custom panels

After library construction, hybridization and capture were performed using the relevant components of Illumina's Nextera Rapid Capture Exome Kit and following the manufacturer's suggested protocol, with the following exceptions: first, all libraries within a library construction plate were pooled prior to hybridization. Second, the Midi plate from Illumina's Nextera Rapid Capture Exome Kit was replaced with a skirted PCR plate to facilitate automation. All hybridization and capture steps were automated on the Agilent Bravo liquid handling system.

Preparation of libraries for cluster amplification and sequencing

After post-capture enrichment, library pools were quantified using qPCR (automated assay on the Agilent Bravo), using a kit purchased from KAPA Biosystems with probes specific to the ends of the adapters. Based on qPCR quantification, pools were normalized to 1.5nM for exome libraries and 1.5nM for genome libraries.

Cluster amplification and sequencing

Cluster amplification of library pools was performed according to the manufacturer's protocol (Illumina) using Exclusion Amplification cluster chemistry and HiSeq X flowcells. Flowcells were sequenced on v2 Sequencing-by-Synthesis chemistry for HiSeq X flowcells. The flowcells are then analyzed using RTA v.2.7.3 or later. Each pool of whole genome libraries was run on paired 151bp runs, reading the dual-indexed sequences to identify molecular indices and sequenced across the number of lanes needed to meet coverage for all libraries in the pool.

Cluster amplification of library pools was performed according to the manufacturer's protocol (Illumina) using Exclusion Amplification cluster chemistry and HiSeq 4000 flowcells. Flowcells were sequenced on v1

Sequencing-by-Synthesis chemistry for HiSeq 4000 flowcells. The flowcells are then analyzed using RTA v.2.7.3 or later. Each pool of whole exome libraries was run on paired 76bp runs, reading the dual-indexed sequences to identify molecular indices and sequenced across the number of lanes needed to meet coverage for all libraries in the pool.

REFERENCES

1. Kapa Biosystems, "KAPA Hyper Prep Kit Technical Data Sheet," KR0961 – v5.16, June 2016.
2. Qiagen, "QIAasympy® DSP Circulating DNA Kit Instructions for Use (Handbook)," Version 1, March 2017.

Exome Express ICE Methods

Library Construction

Library construction was performed as described in Fisher et al., with the following modifications: initial genomic DNA input into shearing was reduced from 3µg to 10-100ng in 50µL of solution. For adapter ligation, Illumina paired end adapters were replaced with palindromic forked adapters, purchased from Integrated DNA Technologies, with unique dual-indexed molecular barcode sequences to facilitate downstream pooling. With the exception of the palindromic forked adapters, the reagents used for end repair, A-base addition, adapter ligation, and library enrichment PCR were purchased from KAPA Biosciences in 96-reaction kits. In addition, during the post-enrichment SPRI cleanup, elution volume was reduced to 30µL to maximize library concentration, and a vortexing step was added to maximize the amount of template eluted.

In-solution hybrid selection

After library construction, hybridization and capture were performed using the relevant components of Illumina's Nextera Rapid Capture Exome Kit and following the manufacturer's suggested protocol, with the following exceptions: first, all libraries within a library construction plate were pooled prior to hybridization. Second, the Midi plate from Illumina's Nextera Rapid Capture Exome Kit was replaced with a skirted PCR plate to facilitate automation. All hybridization and capture steps were automated on the Agilent Bravo liquid handling system.

Preparation of libraries for cluster amplification and sequencing

After post-capture enrichment, library pools were quantified using qPCR (automated assay on the Agilent Bravo), using a kit purchased from KAPA Biosystems with probes specific to the ends of the adapters. Based on qPCR quantification, libraries were normalized to 2nM, then denatured using 0.1 N NaOH on the Hamilton Starlet. After denaturation, libraries were diluted to 20pM using hybridization buffer purchased from Illumina.

Cluster amplification and sequencing

Cluster amplification of denatured templates was performed according to the manufacturer's protocol (Illumina) using HiSeq 4000 cluster chemistry and HiSeq 4000 flowcells. Flowcells were sequenced on v1 Sequencing-by-Synthesis chemistry for HiSeq 4000 flowcells. The flowcells are then analyzed using RTA v.1.18.64 or later. Each pool of whole exome libraries was run on paired 76bp runs, reading the dual-indexed sequences to identify molecular indices and sequenced across the number of lanes needed to meet

coverage for all libraries in the pool.

REFERENCES

1. Fisher S, Barry A, Abreu J, Minie B, Nolan J, Delorey TM, Young G, Berlin AM, Blumenstiel B, Cibulskis K, Friedrich D, Johnson R, Juhn F, Reilly B, Shammas R, Stalker J, Sykes SM, Thompson J, Walsh J, Zimmer A,
2. Zwirk Z, Gabriel S, Nicol R, Nusbaum C. A scalable, fully automated process for construction of sequence-ready human exome targeted capture libraries. *Genome Biology* 2011, 12:R1.

Sequencing Data Analysis

Whole exome sequences were captured using Illumina technology and the sequence data processing and analysis was performed using the Picard and Firehose pipelines at the Broad Institute. The Picard pipeline (<http://picard.sourceforge.net>) was used to produce a BAM file with aligned reads. This includes alignment to the GRCh37 human reference sequence using the BWA aligner^[1] and estimation and recalibration of base quality score with the Genome Analysis Toolkit (GATK)^[2]. All sample pairs passed through the Firehose pipeline were subjected to QC testing to test for any tumor/normal and inter-individual contamination as previously described^[3,4]. The MuTect algorithm was used to identify somatic mutations^[4]. To reduce false positive calls, we additionally analyzed reads covering sites of a putative somatic mutation and realigned them with NovoAlign (www.novocraft.com) and performed additional iteration of MuTect inference on newly aligned BAM files. Furthermore, we filtered for false-positive somatic mutation calls using a panel of normals (PoN), generated from the TCGA dataset and oxoG filter^[5].

Small somatic insertions and deletions were detected using Strelka algorithm^[6], MuTect 2 and Snowman (<https://github.com/walaj/svaba>). Insertions and deletions which were called by two out of the three methods mentioned above were used for analysis.

Somatic mutations including single-nucleotide variants, insertions, and deletions were annotated using Oncotator^[7]. The germline somatic variants were analysed using the HaplotypeCaller module of GATK^[2]. Significance of identified somatic mutations was analyzed using MutSig2CV^[8], which uses patient and gene-specific mutation rates to estimate a background model of predicted mutation incidence across the genome. MutSig2CV then factors in biological co-variables such as replication timing and gene-expression level on a gene-by-gene basis to account for the increased mutational rate of certain classes of genes. To analyze somatic copy number alterations (SCNA) from whole exome data, we used ReCapSeg, which assesses homolog-specific copy ratios from segmental estimates of multipoint allelic copy ratios at heterozygous loci incorporating the statistical phasing software (BEAGLE) and population haplotype panels (HAPMAP3)^[9,10].

For copy number alteration significance analysis, segmented copy number data was analyzed by GISTIC 2.0, to identify significantly recurring focal and arm-level amplification/deletion peaks.^[11] Allele-specific SCNAs and tumor ploidy/purity status were assessed using ABSOLUTE.^[12]

Angiosarcoma tumor samples which had a purity of 10% or more were submitted to cBioPortal (<http://www.cbioportal.org/index.do>).

See Appendix 2, pg. 9, for detailed pipeline information.

Assessment of Tumor Mutation Burden (TMB)

TMB (mutation per megabase) was calculated as the total number of mutations (non-synonymous + synonymous) detected for a given sample divided by the length of the total genomic target region captured with the whole exome sequencing^[13]. Samples with TMB ≥ 10 mutations per megabase were classified as hypermutated.

REFERENCES

1. Li, H. & Durbin, R. Fast and accurate short read alignment with Burrows-Wheeler transform. *Bioinformatics* **25**, 1754–60 (2009).
2. McKenna, A. *et al.* The Genome Analysis Toolkit: a MapReduce framework for analyzing next-generation DNA sequencing data. *Genome Res.* **20**, 1297–303 (2010).
3. Berger, M. F. *et al.* The genomic complexity of primary human prostate cancer. *Nature* **470**, 214–20 (2011).
4. Cibulskis, K. *et al.* Sensitive detection of somatic point mutations in impure and heterogeneous cancer samples. *Nat. Biotechnol.* **31**, 213–9 (2013).
5. Costello, M. *et al.* Discovery and characterization of artifactual mutations in deep coverage targeted capture sequencing data due to oxidative DNA damage during sample preparation. *Nucleic Acids Res.* **41**, e67– (2013).
6. Saunders, C. T. *et al.* Strelka: accurate somatic small-variant calling from sequenced tumor-normal sample pairs. *Bioinformatics* **28**, 1811–7 (2012).
7. Ramos, A. H. *et al.* Oncotator: cancer variant annotation tool. *Hum. Mutat.* **36**, E2423-9 (2015).
8. Lawrence, M. *et al.* Mutational heterogeneity in cancer and the search for new cancer-associated genes. *Nature* **499**, 214-218 (2013)
9. Browning, B. L. & Yu, Z. Simultaneous genotype calling and haplotype phasing improves genotype accuracy and reduces false-positive associations for genome-wide association studies. *Am. J. Hum. Genet.* **85**, 847–861 (2009).
10. International HapMap, C. *et al.* Integrating common and rare genetic variation in diverse human populations. *Nature* **467**, 52–58 (2010).
11. Beroukhi, R. *et al.* Assessing the significance of chromosomal aberrations in cancer: methodology and application to glioma. *Proc. Natl Acad. Sci. USA* **104**, 20007–20012(2007)
12. Carter, S.L. *et al.* Absolute quantification of somatic DNA alterations in human cancer. *Nat. Biotechnol.* **30**, 413–421 (2012).
13. Chalmers *et al.* Analysis of 100,000 human cancer genomes reveals the landscape of tumor mutational burden *Genome Medicine* **9**:34 (2017)

Appendix 1: Data Dictionary

Medical Record Data	
Date of Diagnosis	Date of biopsy with confirmed angiosarcoma
Diagnostic biopsy location	Location in the body where diagnostic biopsy was performed
Primary site	Site of origin of disease
Primary site notes	Notes regarding site of origin of disease
Sex	Male Female
Date of Birth	Date of birth

Pathology Data	
Sample Collection Time	Days from primary diagnosis
Biopsy location	Location in the body where procedure was performed
Biopsy procedure type	Core Biopsy Punch Biopsy Blood Biopsy Mastectomy Resection Excision Other
Vasoformative Characteristics	Yes No
Epithelioid Characteristics	Yes No
Spindle Cell Characteristics	Yes No
Nuclear Grade	High Low

Appendix 2: Commands used to execute the methods

Key	Workflow/Task	Commands
Contest	ContEst for Capture Array-Free	<pre> java -Djava.io.tmpdir=/fh/subscription-XRC/ContaminationAnalysis/ASC-ASCProject_Normal_Sample/27115945/tmp -Xmx512m -jar /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ContaminationAnalysis/broadinstitute.org/cancer.genome.analysis/0026 2/107/Queue-1.4-437-g6b8a9e1-svn-35362.jar -S /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ContaminationAnalysis/broadinstitute.org/cancer.genome.analysis/0026 2/107/ContaminationPipeline.scala -reference /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta -interval /xchip/cga/reference/hg19/whole_exome_agilent_1.1_refseq_plus_3_boosters_plus_10bp_padding_minus_mito.Homo_sapiens_assembly19.targets.in terval_list-out ASC-ASCProject_Normal_Sample.contamination.txt-bam/seq/picard_aggregation/RP- 1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam -nbam /seq/picard_aggregation/RP- 1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam-array none-pop /xchip/cga/reference/hg19/hg19_population_stratified_af_hapmap_3.3.fixed.vcf-faf true -run -array_interval /xchip/cga/reference/hg19/SNP6.hg19.interval_list </pre>
SNV	Run Mutect with Realignment Filter Pairs NovoAlign and Agilent & ICE (Capture-Pair) for FFPE or FF without ExAC CSQ	<p>Call Somatic Mutations for Capture Workflow without ExAC CSQ</p> <p>1.Call Somatic Mutations for Capture (configure) -</p> <pre> sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallSomaticMutations/broadinstitute.org/cancer.genome.analysis/00004 /131/runBroadJava7.sh java -Xmx2g -jar /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallSomaticMutations/broadinstitute.org/cancer.genome.analysis/00004 /131/muTect-1.1.6.jar --read_group_black_list /xchip/cga/reference/read_group_blacklist.v1.txt --analysis_type MuTect --intervals /fh/subscription- XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/27003052/iteration1/gatk-scatter.0000000043.interval_list --normal_sample_name ASC-ASCProject_Normal_Sample-I.normal/seq/picard_aggregation/RP- 1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam --tumor_sample_name ASC-ASCProject_Tumor_Sample-I.tumor /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam --reference_sequence /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta --dbsnp /xchip/cga/reference/hg19/dbsnp_134_b37.leftAligned.vcf -- cosmic /xchip/cga/reference/hg19/hg19_cosmic_v54_120711.vcf --normal_panel /xchip/cga/reference/hg19/refseq_exome_10bp_hg19_300_1kg_normal_panel.vcf --out ASC-ASCProject_Tumor_Normal.call_stats.txt -- coverage_file ASC-ASCProject_Tumor_Normal.coverage.wig.txt --power_file ASC-ASCProject_Tumor_Normal.power.wig.txt -- downsample_to_coverage 10000 --enable_extended_output --fraction_contamination 0.007 </pre> <p>2.CallStats to MAFLite for Capture (configure) -</p> <pre> /usr/bin/perl /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallStatsToMaflite/broadinstitute.org/cancer.genome.analysis/00162/14/ call_stats_to_maflite.pl /fh/subscription-XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/27003052/iteration1/ASC- ASCProject_Tumor_Normal.call_stats.txt 37FSTARASC-ASCProject_Tumor_Normal.maf tumor_f_init_t_lod_t_lod_fstar_t_alt_count_t_ref_count_judgement </pre> <p>3.Summarize Somatic Coverage for Capture (configure) -</p> <pre> perl /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/SummarizeWigFile/broadinstitute.org/cancer.genome.analysis/10246/1/ summarizeWigFile.pl /fh/subscription-XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/27003052/iteration1/ASC- ASCProject_Tumor_Normal.coverage.wig.txt.gz ASC-ASCProject_Tumor_Normal.somatic_coverage_summary.txt </pre> <p>Oncotator Capture SNP Workflow without ExAC CSQ</p> <p>4.Apply SNP Maflite Validation for Capture (configure)</p> <pre> /broad/software/free/Linux/redhat_5_x86_64/pkgs/sun-java-jdk_1.6.0-21_x86_64/bin/java -Xmx1g -jar /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ApplyMAFValidation/broadinstitute.org/cancer.genome.analysis/00163/ 23/ApplyMAFValidation.jar M=/fh/subscription-XRC/CallStatsToMaflite/ASC-ASCProject_Tumor_Normal/27022277/ASC- ASCProject_Tumor_Normal.mafOUTPUT_MAF=ASC-ASCProject_Tumor_Sample.maf.annotatedMATCH_MODE=SampleV=/cga/tcga- gsc/svnreference/validation </pre> <p>5.Oncotator SNP for Capture (configure)</p> <pre> sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Oncotator_v1/broadinstitute.org/cancer.genome.analysis/10202/68/onco tator.sh MAFLITE TCGAMAF /fh/subscription-XRC/ApplyMAFValidation/ASC-ASCProject_Tumor_Normal/27022407/ASC- ASCProject_Tumor_Sample.maf.annotated ASC-ASCProject_Tumor_Normal.snp.capture.maf.annotated hg19 /xchip/cga/reference/annotation/db/oncotator_v1_ds_Sept292015/ /xchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Oncotator_v1/broadinstitute.org/cancer.genome.analysis/10202/68/ EFFECT /xchip/tcga/Tools/python_fh/python_fh_env_March302015 --log_name oncotator_firehose.log --prepend </pre> <p>6.Oncotator - Drop ExAC CSQ field (configure)</p> <pre> python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/SubsetMafColumns/broadinstitute.org/cancer.genome.analysis/10917/1 6/subset_maf_by_columns.py --maf /fh/subscription-XRC/Oncotator_v1/ASC-ASCProject_Tumor_Normal/27022437/ASC- ASCProject_Tumor_Normal.snp.capture.maf.annotated --columns /home/unix/dkim/reference/firehose/oncotator_whitelist_no_ExAC_CSQ.txt --out ASC-ASCProject_Tumor_Normal.snp.capture.maf.annotated.no_exac_csq </pre> <p>7.PoN filter (configure)</p> <pre> /broad/software/free/Linux/redhat_5_x86_64/pkgs/python_2.7.1-sqlite3-rtrees/bin/python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/PoN_filter_py/broadinstitute.org/cancer.genome.analysis/11041/16/PO N_filter_fh.py -m /fh/subscription-XRC/Oncotator_v1/ASC-ASCProject_Tumor_Normal/30662860/ASC- ASCProject_Tumor_Normal.snp.capture.maf.annotated -p /cga/fh/pancan_data/pon/pon12/final_summed_tokens.hist.bin -g 0.005 -t 5 -wc 0.5 -pc 0.5 -a 0.2 -r TRUE -o ASC-ASCProject_Tumor_Normal.PoNfiltered.maf-n l -l -2.5 -w c 0.5 --use_static </pre> <p>Capture OxoG Filter v3 Workflow without ExAC CSQ</p> <p>8.annotate Picard OxoQ for Capture (configure)</p> <pre> sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/annotatePicardOxoQ/broadinstitute.org/cancer.genome.analysis/10360/6 /annotatePicardOxoQ.sh -i ASC-ASCProject_Tumor_Normal-b /seq/picard_aggregation/RP- 1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam -r /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta -d /xchip/cga/reference/hg19/dbsnp_134_b37.leftAligned.vcf-c CCG-o. </pre> <p>9.Create OxoG Interval List for Capture (configure)</p> <pre> sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/createOxoGIntervalList/broadinstitute.org/cancer.genome.analysis/1006 7/22/createOxoGIntervals.sh /fh/subscription-XRC/Oncotator_v1/ASC-ASCProject_Tumor_Normal/27022437/ASC- </pre>

```

ASCProject_Tumor_Normal.snp.capture.maf.annotated ASC-ASCProject_Tumor_Normal.oxoG.interval_list
10. Append Picard OxoQ to Maf for Capture without ExAC CSQ (configure)
sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/AppendAnnotation2MAF/broadinstitute.org/cancer.genome.analysis/10
363/16/AppendAnnotation2MAF.sh -i ASC-ASCProject_Tumor_Normal -m /fh/subscription-XRC/SubsetMafColumns/ASC-
ASCProject_Tumor_Normal/27023157/ASC-ASCProject_Tumor_Normal.snp.capture.maf.annotated.no_exac_csq_lite.maf -f picard_oxoQ -v 52.52 -
o
11. Create OxoG Metrics for Capture (configure)
java -Xmx2g -jar
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/oxoGMetrics/broadinstitute.org/cancer.genome.analysis/10069/8/Geno
meAnalysisTK.jar --analysis_type OxoGMetrics -R /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta -I
/seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam -L /fh/subscription-
XRC/createOxoGIntervalList/ASC-ASCProject_Tumor_Normal/27072480/ASC-ASCProject_Tumor_Normal.oxoG.interval_list -o ASC-
ASCProject_Tumor_Normal.oxoG.metrics.txt
12. Append OxoG Information to MAF for Capture (configure)
sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/appendOxoGInfo/broadinstitute.org/cancer.genome.analysis/10138/8/a
ppendOxoGInfo.sh --onlyAddColumnsToCopy /fh/subscription-XRC/oxoGMetrics/ASC-ASCProject_Tumor_Normal/27072528/ASC-
ASCProject_Tumor_Normal.oxoG.metrics.txt /fh/subscription-XRC/AppendAnnotation2MAF/ASC-ASCProject_Tumor_Normal/27072507/ASC-
ASCProject_Tumor_Normal.picard_oxoQ.maf.annotated ASC-ASCProject_Tumor_Normal.oxoGInfo.maf.annotated
13. Filter OxoG Artifact Third Incarnation for Capture (configure)
sh /xchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2013a
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/oxoGFilter_v3/broadinstitute.org/cancer.genome.analysis/10139/69/ --
with-display
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/oxoGFilter_v3/broadinstitute.org/cancer.genome.analysis/10139/69/start
FilterMAFFile /fh/subscription-XRC/appendOxoGInfo/ASC-ASCProject_Tumor_Normal/27115529/ASC-
ASCProject_Tumor_Normal.oxoGInfo.maf.annotated ASC-ASCProject_Tumor_Normal.oxoG3.maf.annotated / 0 1 0.96 0.01 -1 36 1.5
14. Generate OxoGv3 Report for Capture (configure)
Rscript
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/oxoGFilter_v3Report/broadinstitute.org/cancer.genome.analysis/10153/
41/createReport.R -L /xchip/tcga/Tools/Nozzle/v1.current -C /fh/subscription-XRC/oxoGFilter_v3/ASC-
ASCProject_Tumor_Normal/27115552/allCases.txt -F /fh/subscription-XRC/oxoGFilter_v3/ASC-ASCProject_Tumor_Normal/27115552/figures -
T /fh/subscription-XRC/oxoGFilter_v3/ASC-ASCProject_Tumor_Normal/27115552/caseTableData.tsv -M /fh/subscription-XRC/oxoGFilter_v3/ASC-
ASCProject_Tumor_Normal/27115552/ASC-ASCProject_Tumor_Normal.oxoG3.maf.annotated -R /fh/subscription-XRC/oxoGFilter_v3/ASC-
ASCProject_Tumor_Normal/27115552/ASC-ASCProject_Tumor_Normal.oxoG3.maf.annotated.all.maf.annotated

Capture Filter FFPE Workflow
15. Annotate FFPE bias for Capture (configure)
sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/annotateFFPEbias/broadinstitute.org/cancer.genome.analysis/10653/7/a
nnotatFFPEbias.sh -i ASC-ASCProject_Tumor_Normal -b /seq/picard_aggregation/RP-
1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam -r /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta
-d /xchip/cga/reference/hg19/dbsnp_134_b37.leftAligned.vcf -f G -t A -c CG -x ffpe_metrics -s 40000000 -o
16. Append ffpeQ to Maf for Capture (configure)
sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/AppendAnnotation2MAF/broadinstitute.org/cancer.genome.analysis/10
363/16/AppendAnnotation2MAF.sh -i ASC-ASCProject_Tumor_Normal -m /fh/subscription-XRC/oxoGFilter_v3/ASC-
ASCProject_Tumor_Normal/27115552/ASC-ASCProject_Tumor_Normal.oxoG3.maf.annotated -f ffpe_Q -v 31.50 -o
17. Filter FFPE orientation bias for capture (configure)
sh /xchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2012a
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/orientationBiasFilter/broadinstitute.org/cancer.genome.analysis/10884/1
3/ --with-display
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/orientationBiasFilter/broadinstitute.org/cancer.genome.analysis/10884/1
3/orientationBiasFilter /fh/subscription-XRC/AppendAnnotation2MAF/ASC-ASCProject_Tumor_Normal/27115594/ASC-
ASCProject_Tumor_Normal.ffpe_Q.maf.annotated ASC-ASCProject_Tumor_Normal.ffpeBias.maf.annotated / 0 1 0.96 0.01 -1 30 1.5 G A i_ffpe
18. ffpe filter Report (configure)
Rscript
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ffpeReport/broadinstitute.org/cancer.genome.analysis/10661/6/ffpeRep
ort.R -L /xchip/tcga/Tools/Nozzle/v1.current -C /fh/subscription-XRC/orientationBiasFilter/ASC-ASCProject_Tumor_Normal/27115864/allCases.txt -
F /fh/subscription-XRC/orientationBiasFilter/ASC-ASCProject_Tumor_Normal/27115864/figures/ -T /fh/subscription-XRC/orientationBiasFilter/ASC-
ASCProject_Tumor_Normal/27115864/caseTableData.tsv -M /fh/subscription-XRC/orientationBiasFilter/ASC-
ASCProject_Tumor_Normal/27115864/ASC-ASCProject_Tumor_Normal.ffpeBias.maf.annotated -R /fh/subscription-XRC/orientationBiasFilter/ASC-
ASCProject_Tumor_Normal/27115864/ASC-ASCProject_Tumor_Normal.ffpeBias.unfiltered.maf.annotated

Call Somatic Mutations for Capture Workflow Realignment Filter Pairs NovoAlign
19. Realignment Filter Paired Reads novoalign (FFPE or FF)
sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReadAnalysis_QC_PairedReads_bwamem/broadinstitute.org/cancer.gen
ome.analysis/10845/21/readalignment_wrapper.sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReadAnalysis_QC_PairedReads_bwamem/broadinstitute.org/cancer.gen
ome.analysis/10845/21/ /fh/subscription-XRC/orientationBiasFilter/ASC-ASCProject_Tumor_Normal/27115864/ASC-
ASCProject_Tumor_Normal.ffpeBias.maf.annotated /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam
/seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.SALIVA.bam
/seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta
/xchip/cga_home/mara/projects/readalignment/novoalign/hg19.decoy.nix ASC-ASCProject_Tumor_Sample ASC-ASCProject_Normal_Sample
ASC-ASCProject_Tumor_Normal realigned bam mutation interval_list length_wex_pairs_novo
bamfile_mutect_realignedreads_control_wex_pairs_novobamfile_mutect_realignedreads_case_wex_pairs_novo novoaligndebug
20. Call Somatic Mutations for WEX Pairs Novoalign Realign
sh
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallSomaticMutations/broadinstitute.org/cancer.genome.analysis/00004
/131/runBroadJava7.sh java -Xmx2g -jar
/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallSomaticMutations/broadinstitute.org/cancer.genome.analysis/00004
/131/muTect-1.1.6.jar --analysis_type MuTect --intervals /fh/subscription-XRC/CallSomaticMutations/ASC-
ASCProject_Tumor_Normal/30661616/iteration1/gatk-scatter.000000009.interval_list --normal_sample_name ASC-ASCProject_Normal_Sample -
Inormal /fh/subscription-XRC/ReadAnalysis_QC_PairedReads_bwamem/ASC-ASCProject_Tumor_Normal/30657886/ASC-
ASCProject_Normal_Sample-final_filtered_control.bam --tumor_sample_name ASC-ASCProject_Tumor_Sample-I.tumor /fh/subscription-
XRC/ReadAnalysis_QC_PairedReads_bwamem/ASC-ASCProject_Tumor_Normal/30657886/ASC-ASCProject_Tumor_Sample-
final_filtered_case.bam --reference_sequence /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta --dbsnp
/xchip/cga/reference/hg19/dbsnp_134_b37.leftAligned.vcf --cosmic /xchip/cga/reference/hg19/hg19_cosmic_v54_120711.vcf --out ASC-
ASCProject_Tumor_Normal.call_stats.txt --coverage_file ASC-ASCProject_Tumor_Normal.coverage.wig.txt --power_file ASC-
ASCProject_Tumor_Normal.power.wig.txt --downsample_to_coverage 9999999 --enable_extended_output --fraction_contamination 0.007 --
force_output

```

		<p>21. Realignment Filter Call Stats WEX Pairs novo</p> <pre>sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReadAnalysis_QC_CallStatsFilter/broadinstitute.org/cancer.genome.analysis/10215/9/callstats_filter_wrapper.sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReadAnalysis_QC_CallStatsFilter/broadinstitute.org/cancer.genome.analysis/10215/9/ /fh/subscription-XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/27003052/iteration1/ASC-ASCProject_Tumor_Normal.call_stats.txt /fh/subscription-XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/30661616/iteration1/ASC-ASCProject_Tumor_Normal.call_stats.txt ASC-ASCProject_Tumor_Normal</pre> <p>22. CallStats to MAF Lite for Capture Realign Novo Pairs</p> <pre>/usr/bin/perl /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallstatsToMafLite/broadinstitute.org/cancer.genome.analysis/00162/14/call_stats_to_mafLite.pl /fh/subscription-XRC/ReadAnalysis_QC_CallStatsFilter/ASC-ASCProject_Tumor_Normal/30662397/ASC-ASCProject_Tumor_Normal.filtered.annotated.callstats.txt 37 FSTAR ASC-ASCProject_Tumor_Normal.maf tumor_f_init t lod t lod fstar t alt count t ref count judgement</pre> <p>23. Apply SNP MAF Lite Validation for Capture Realign Pairs Novo</p> <pre>/broad/software/free/Linux/redhat_5_x86_64/pkgs/sun-java-jdk_1.6.0-21_x86_64/bin/java -Xmx1g -jar /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ApplyMAFValidation/broadinstitute.org/cancer.genome.analysis/00163/23/ApplyMAFValidation.jar M=/fh/subscription-XRC/CallstatsToMafLite/ASC-ASCProject_Tumor_Normal/30662533/ASC-ASCProject_Tumor_Normal.maf OUTPUT_MAF=ASC-ASCProject_Tumor_Sample.maf.annotated MATCH_MODE=Sample V=/cga/tcga-gsc/svreference/validation</pre> <p>24. Oncotator SNP for Capture Realign Pairs Novo</p> <pre>sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Oncotator_v1/broadinstitute.org/cancer.genome.analysis/10202/70/oncotator.sh MAF_LITE_TCGAMAF /fh/subscription-XRC/ApplyMAFValidation/ASC-ASCProject_Tumor_Normal/30662698/ASC-ASCProject_Tumor_Sample.maf.annotated ASC-ASCProject_Tumor_Normal.snp.capture.maf.annotated hg19 /xcchip/cga/reference/annotation/db/oncotator_v1_ds_Sept292015/ /xcchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Oncotator_v1/broadinstitute.org/cancer.genome.analysis/10202/70/EFFECT /xcchip/tcga/Tools/python_fh/python_fh_env_March302015 --log_name oncotator_firehose.log --prepend</pre> <p>25. PoN Filter (Agilent) - Capture Realignment Filter Pairs NovoAlign</p> <pre>/broad/software/free/Linux/redhat_5_x86_64/pkgs/python_2.7.1-sqlite3-rtrees/bin/python /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/PoN_filter_py/broadinstitute.org/cancer.genome.analysis/11041/22/PO_N_filter_fh.py -m /fh/subscription-XRC/Oncotator_v1/ASC-ASCProject_Tumor_Normal/30662860/ASC-ASCProject_Tumor_Normal.snp.capture.maf.annotated -p /cga/fh/pancan_data/pon/pon9/final_summed_tokens.hist.bin -g 0.005 -t 5 -wc 0.5 -pc 0.5 -a 0.2 -r TRUE -o ASC-ASCProject_Tumor_Normal.novoalign.PoNfiltered_agilent.maf -n 1 -l -2.5 -w_c 0.5 --use_static</pre> <p>26. PoN Filter (ICE from Agilent) - Capture Realignment Filter Pairs NovoAlign</p> <pre>/broad/software/free/Linux/redhat_5_x86_64/pkgs/python_2.7.1-sqlite3-rtrees/bin/python /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/PoN_filter_py/broadinstitute.org/cancer.genome.analysis/11041/24/PO_N_filter_fh.py -m /fh/subscription-XRC/PoN_filter_py/ASC-ASCProject_Tumor_Normal/30688082/ASC-ASCProject_Tumor_Normal.novoalign.PoNfiltered_agilent.maf -p /cga/fh/pancan_data/pon/pon12/final_summed_tokens.hist.bin -g 0.005 -t 5 -wc 0.5 -pc 0.5 -a 0.2 -r TRUE -o ASC-ASCProject_Tumor_Normal.novoalign.PoNfiltered_agilent_ICE.maf -n 1 -l -2.5 -w_c 0.5 --use_static --no_pon_columns</pre>
Indel (1)	Run Strelka with MafMaster Filter (Capture-Pair)	<p>1. Run Strelka on Capture for Pairs -</p> <pre>sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/RunStrelka/broadinstitute.org/cancer.genome.analysis/10439/11/runStrelka.sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/RunStrelka/broadinstitute.org/cancer.genome.analysis/10439/11/seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta ASC-ASCProject_Tumor_Normal strelka_config_bwa_cga_exome.ini</pre> <p>2. Annotate Pass Strelka Indel VCF for Pairs -</p> <pre>python /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorForStrelka/broadinstitute.org/cancer.genome.analysis/10543/8/oncotator/Oncotator.py -v --input_format=VCF --output_format=TCGAMAF --db-dir=/xcchip/cga/reference/annotation/db/oncotator_v1_ds_April052016 --no-multicore --default_config /xcchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /fh/subscription-XRC/FixStrelkaVCFForAnnotation/ASC-ASCProject_Tumor_Normal/27011318/tumorOnly.fixed.vcf ASC-ASCProject_Tumor_Normal.strelka.pass.somatic.indels.maf hg19-a center:broad.mit.edu -a individual_barcode:ASC-ASCProject_pt_id -d normal_barcode:ASC-ASCProject_Normal_Sample -a platform:illumina -a normal_uuid:NA -a tumor_barcode:ASC-ASCProject_Tumor_Sample -a tumor_uuid:NA -a tumor_subtype:Metastatic</pre> <p>3. Annotate Pass Strelka SNV VCF for Pairs -</p> <pre>python /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorForStrelka/broadinstitute.org/cancer.genome.analysis/10543/8/oncotator/Oncotator.py -v --input_format=VCF --output_format=TCGAMAF --db-dir=/xcchip/cga/reference/annotation/db/oncotator_v1_ds_April052016 --no-multicore --default_config /xcchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /fh/subscription-XRC/FixStrelkaVCFForAnnotation/ASC-ASCProject_Tumor_Normal/27011317/tumorOnly.fixed.vcf ASC-ASCProject_Tumor_Normal.strelka.pass.somatic.snvs.maf hg19-a center:broad.mit.edu -a individual_barcode:ASC-ASCProject_pt_id -d normal_barcode:ASC-ASCProject_Normal_Sample -a platform:illumina -a normal_uuid:NA -a tumor_barcode:ASC-ASCProject_Tumor_Sample -a tumor_uuid:NA -a tumor_subtype:Metastatic</pre> <p>4. Strelka SNV fix INDEL allele counts -</p> <pre>sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/strelkamaf2tcgamaf/broadinstitute.org/cancer.genome.analysis/11325/1/strelkamaf2tcgamaf.sh -i ASC-ASCProject_Tumor_Normal-m /fh/subscription-XRC/AnnotateStrelkaVCF/ASC-ASCProject_Tumor_Normal/31666346/OncotatorForStrelka_31666749/ASC-ASCProject_Tumor_Normal.strelka.pass.somatic.indels.maf -x strelka.fix.indel.maf -o .</pre> <p>5. Strelka SNV fix SNV allele counts -</p> <pre>sh /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/strelkamaf2tcgamaf/broadinstitute.org/cancer.genome.analysis/11325/1/strelkamaf2tcgamaf.sh -i ASC-ASCProject_Tumor_Normal-m /fh/subscription-XRC/AnnotateStrelkaVCF/ASC-ASCProject_Tumor_Normal/31666342/OncotatorForStrelka_31666768/ASC-ASCProject_Tumor_Normal.strelka.pass.somatic.snvs.maf -x strelka.fix.snv.maf -o .</pre> <p>6. Combine SNP and Indel Strelka MAF for Capture -</p> <pre>python /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/tsvCatFiles/broadinstitute.org/cancer.genome.analysis/10155/3/tsvConcatFiles.py /fh/subscription-XRC/strelkamaf2tcgamaf/ASC-ASCProject_Tumor_Normal/31667385/output/ASC-ASCProject_Tumor_Normal.strelka.fix.snv.maf /fh/subscription-XRC/strelkamaf2tcgamaf/ASC-ASCProject_Tumor_Normal/31667383/output/ASC-ASCProject_Tumor_Normal.strelka.fix.indel.maf --outputFilename=ASC-ASCProject_Tumor_Normal.combined.maf</pre> <p>7. PoN filter indels -</p> <pre>sh /xcchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2012a /xcchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MafMasterFilter/broadinstitute.org/cancer.genome.analysis/10822/1/</pre>

		<pre> /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MafMasterFilter/broadinstitute.org/cancer.genome.analysis/10822/1/master_filter_wrapper /fh/subscription-XRC/AnnotateStrelkaVCF/ASC-ASCProject_Tumor_Normal/31666346/OncotatorForStrelka_31666749/ASC-ASCProject_Tumor_Normal.strelka.pass.somatic.indels.maf /cga/fh/pancan_data/pon/pon10/final_summed_tokens.hist.bin ASC-ASCProject_Tumor_Normal 8.Strelka PoN Filter (ICE from Agilent) - sh /xchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2012a /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MafMasterFilter/broadinstitute.org/cancer.genome.analysis/10822/1/ /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MafMasterFilter/broadinstitute.org/cancer.genome.analysis/10822/1/master_filter_wrapper /fh/subscription-XRC/MafMasterFilter/ASC-ASCProject_Tumor_Normal/31668062/ASC-ASCProject_Tumor_Normal.pon_filtered.txt /cga/fh/pancan_data/pon/pon12/final_summed_tokens.hist.bin ASC-ASCProject_Tumor_Normal 10.Aggregate Strelka Indels -ICE/Agilent filtered - python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/tsvCat/broadinstitute.org/cancer.genome.analysis/10083/10/tsvConcatLstFile.py /fh/subscription-XRC/tsvCat/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33627161/inputListFile.input.tsv PR_Wagle_ASC_Pairs_Capture_All_Pairs.strelka_indels_after_ccpm_ice_agilent </pre>
Indel (2)	M2	<pre> 1. M2 Scatter - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/M2/broadinstitute.org/cancer.genome.analysis/10855/28/runBroadJava7.sh java -Xmx4g -jar /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/M2/broadinstitute.org/cancer.genome.analysis/10855/28/GenomeAnalysisTK_latest_unstable.jar --analysis_type M2 --intervals /fh/subscription-XRC/M2/ASC-ASCProject_Tumor_Normal/32671507/iteration1/./gatk-scatter.0000000020.interval_list -l normal /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam -l tumor/seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam --reference_sequence /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta --contamination 0.007 --dbsnp /xchip/cga/home/kcibul/analysis/dream/mutect_ref/dbsnp_138.b37.vcf --cosmic /xchip/cga/reference/hg19/hg19_cosmic_v54_120711.vcf --normal_panel /cga/tcga-gsc/mutect/panel_of_normals/panel_of_normals_m2_paad/agilent hg19_m2_paad_149_normal_panel.vcf --out ASC-ASCProject_Tumor_Normal.full.vcf 2.VCF Filter to Pass - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/VCF_filterpass/broadinstitute.org/cancer.genome.analysis/10856/8/wrapper.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/VCF_filterpass/broadinstitute.org/cancer.genome.analysis/10856/8/ /fh/subscription-XRC/M2/ASC-ASCProject_Tumor_Normal/32671507/ASC-ASCProject_Tumor_Normal.full.vcf ASC-ASCProject_Tumor_Normal /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta 3.Oncotator for M2 - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorForStrelka/broadinstitute.org/cancer.genome.analysis/10543/8 /oncotator.sh VCF TCGAMAF /fh/subscription-XRC/VCF_filterpass/ASC-ASCProject_Tumor_Normal/32684786/ASC-ASCProject_Tumor_Normal.pass.vcf ASC-ASCProject_Tumor_Normal.m2_pass.maf hg19 /xchip/cga/reference/annotation/db/oncotator_v1_ds_Sep1292015/ /xchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorForStrelka/broadinstitute.org/cancer.genome.analysis/10543/8 /center:broad.mit.edu/individual_barcode:ASC-ASCProject_pt_idnormal_barcode:ASC-ASCProject_Normal_Sampleplatform:illumina normal_uuid:ASC-ASCProject_Normal_Sample tumor_barcode:ASC-ASCProject_Tumor_Sampletumor_uuid:ASC-ASCProject_Tumor_Sampletumor_subtype:NA /xchip/tcga/Tools/python fh/python fh env_March302015/ 4.Add Allelic Fraction Columns to MAF - python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/M2_reformat_maf/broadinstitute.org/cancer.genome.analysis/10997/1/ make_oxoG_compatible.py --maf /fh/subscription-XRC/OncotatorForStrelka/ASC-ASCProject_Tumor_Normal/32684814/ASC-ASCProject_Tumor_Normal.m2_pass.maf --out ASC-ASCProject_Tumor_Normal </pre>
Indel (3)	Snowman	<pre> 1.SnowmanWithIndel_Capture - /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Snowman/broadinstitute.org/cancer.genome.analysis/10982/134/snow.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Snowman/broadinstitute.org/cancer.genome.analysis/10982/134/ /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam ASC-ASCProject_Tumor_Normal 0 -n /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam -p4 -D/xchip/gistic/Jeremiah/Projects/SnowmanFilters/dbsnp_138.b37_indel.vcf -B/xchip/gistic/Jeremiah/Projects/SnowmanFilters/HengLiMask/snowman_blacklist.bed -Y/xchip/gistic/Jeremiah/Projects/SnowmanFilters/viral.1.1.genomic_ns.fna 2.Oncotator Snowman Indel for Capture - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorVCFGZ/broadinstitute.org/cancer.genome.analysis/10977/4/oncotator.sh VCF TCGAMAF /fh/subscription-XRC/Snowman/ASC-ASCProject_Tumor_Normal/31865571/ASC-ASCProject_Tumor_Normal.snowman.somatic.indel.vcf ASC-ASCProject_Tumor_Normal.indel.capture.maf.annotated hg19 /xchip/cga/reference/annotation/db/oncotator_v1_ds /xchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorVCFGZ/broadinstitute.org/cancer.genome.analysis/10977/4/ CANONICAL /xchip/tcga/Tools/python fh/python fh env_Dec122014/ --log_name oncotator_firehose.log --prepend --infer-onps -c /xchip/cga/reference/annotation/db/tx_exact_uniprot_matches.txt 3.Oncotator Snowman Indel for Capture unfiltered - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorVCFGZ/broadinstitute.org/cancer.genome.analysis/10977/4/oncotator.sh VCF TCGAMAF /fh/subscription-XRC/Snowman/ASC-ASCProject_Tumor_Normal/31865571/ASC-ASCProject_Tumor_Normal.snowman.unfiltered.somatic.indel.vcf ASC-ASCProject_Tumor_Normal.indel.capture.unfiltered.maf.annotated hg19 /xchip/cga/reference/annotation/db/oncotator_v1_ds /xchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/OncotatorVCFGZ/broadinstitute.org/cancer.genome.analysis/10977/4/ CANONICAL /xchip/tcga/Tools/python fh/python fh env_Dec122014/ --log_name oncotator_firehose.log --prepend --infer-onps -c /xchip/cga/reference/annotation/db/tx_exact_uniprot_matches.txt 4.Snowman Annotate Capture - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/SnowmanAnnotate/broadinstitute.org/cancer.genome.analysis/11109/25 /run.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/SnowmanAnnotate/broadinstitute.org/cancer.genome.analysis/11109/25 /snow-annotate.R -g TRUE -W 5 -H 5 -s 0 -i /fh/subscription-XRC/Snowman/ASC-ASCProject_Tumor_Normal/31865571/ASC-ASCProject_Tumor_Normal.snowman.somatic.vcf -o ASC-ASCProject_Tumor_Normal -d 8 -t 4 5.Aggregate Individual MAF Indel Snowman - python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/tsvCat/broadinstitute.org/cancer.genome.analysis/10083/10/tsvConcatL </pre>

		<pre> istFile.py /fh/subscription-XRC/tsvCat/PR_Wagle_ASC_Pairs_Capture_All_Pairs/32757023/input.ListFile.input.tsv PR_Wagle_ASC_Pairs_Capture_All_Pairs.indel.snowman.maf.annotated 6.Aggregate Individual MAF Indel Snowman unfiltered - python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/tsvCat/broadinstitute.org/cancer.genome.analysis/10083/10/tsvConcatL istFile.py /fh/subscription-XRC/tsvCat/PR_Wagle_ASC_Pairs_Capture_All_Pairs/32757022/input.ListFile.input.tsv PR_Wagle_ASC_Pairs_Capture_All_Pairs.indel.snowman.unfiltered.maf.annotated </pre>
ReCapSeq	ReCapSeg_w rokflow	<pre> 1.Picard Target Mapper - python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/PicardTargetMapper/broadinstitute.org/cancer.genome.analysis/11037/ 3/picard_target_mapper.py ASC-ASCProject_Normal_Sample /xchip/cga/reference/seq_target_to_annotations.v4.json /seq/picard_aggregation/RP- 1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam output.txt 2.ReCapSeg Coverage - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReCapSegCoverage/broadinstitute.org/cancer.genome.analysis/10608/ 20/wrapper.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReCapSegCoverage/broadinstitute.org/cancer.genome.analysis/10608/ 20/ -- jar /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ReCapSegCoverage/broadinstitute.org/cancer.genome.analysis/10608/ 20/coverage.jar --sample_id ASC-ASCProject_Normal_Sample --bam /seq/picard_aggregation/RP- 1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam --bed /xchip/cga/reference/hg19/recapseg/target_beds/CRSP_ICE_hg19_wex_illumina_v1.no_X_Y_MT.bed --reference /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta 3.ReCapSeg Read Group Collapse and Proportional Coverage - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/recapseg_pc_collapse/broadinstitute.org/cancer.genome.analysis/11060 /1/recapseg_pc_collapse.sh /fh/subscription-XRC/ReCapSegCoverage/ASC-ASCProject_Normal_Sample/27010468/ASC- ASCProject_Normal_Sample.coverage /xchip/cga/reference/hg19/recapseg/target_beds/CRSP_ICE_hg19_wex_illumina_v1.no_X_Y_MT.bed /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam ASC-ASCProject_Normal_Sample /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/recapseg_pc_collapse/broadinstitute.org/cancer.genome.analysis/11060 /1/ /xchip/tcga/Tools/recapseg/releases/recapseg_env_FH_May182015/ 1.0 4.ReCapSeg Filter TN Segment - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/recapseg_tumor_pcov/broadinstitute.org/cancer.genome.analysis/1066 5/34/recapseg_tumor_pcov.sh /xchip/cga/reference/hg19/recapseg/pon/ICE_Combined_EEW_EWS_v1.pon /fh/subscription- XRC/recapseg_pc_collapse/ASC-ASCProject_Normal_Sample/27010964/ASC-ASCProject_Normal_Sample.pcov /fh/subscription- XRC/recapseg_pc_collapse/ASC-ASCProject_Normal_Sample/27010964/ASC-ASCProject_Normal_Sample.pcov.cr.stat ASC- ASCProject_Normal_Samplerecapseg /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/recapseg_tumor_pcov/broadinstitute.org/cancer.genome.analysis/1066 5/34/ /xchip/tcga/Tools/recapseg/releases/recapseg_env_FH_May182015/ /xchip/cga/reference/hg19/recapseg/target_beds/CRSP_ICE_hg19_wex_illumina_v1.no_X_Y_MT.bed --is-plotting 5.ReCapSeg Caller - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/recapseg_caller/broadinstitute.org/cancer.genome.analysis/10797/10/re capseg_caller.sh /fh/subscription-XRC/recapseg_tumor_pcov/ASC-ASCProject_Normal_Sample/27010966/recapseg.tn.ASC- ASCProject_Normal_Sample.tsv /fh/subscription-XRC/recapseg_tumor_pcov/ASC-ASCProject_Normal_Sample/27010966/recapseg.ASC- ASCProject_Normal_Sample.seg recapseg.ASC-ASCProject_Normal_Sample.called /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/recapseg_caller/broadinstitute.org/cancer.genome.analysis/10797/10/ /xchip/tcga/Tools/recapseg/releases/recapseg_env_FH_May182015/ 6.ReCapSeg Oncotate SCNA Seg File - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Oncotator_v1/broadinstitute.org/cancer.genome.analysis/10202/60/onc otator.sh SEG_FILE.GENE_LIST /fh/subscription-XRC/recapseg_caller/ASC-ASCProject_Normal_Sample/27010973/recapseg.ASC- ASCProject_Normal_Sample.called ASC-ASCProject_Normal_Sample.called.seg.annotated hg19 /xchip/cga/reference/annotation/db/oncotator_v1_ds_June112014/ /xchip/cga/reference/annotation/db/tcgaMAFManualOverrides2.4.config /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Oncotator_v1/broadinstitute.org/cancer.genome.analysis/10202/60/ CANONICAL /xchip/tcga/Tools/python_fh/python_fh_env_June252014/ --log_name oncotator_firehose_seg.log </pre>
Gistic2	Gistic2_Fro m_ReCapSeg _Workflow	<pre> 1.Aggregate_ReCapSeg_Segs - python /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/tsvCat/broadinstitute.org/cancer.genome.analysis/10083/10/tsvConcat ListFile.py /fh/subscription-XRC/tsvCat/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33645478/input.ListFile.input.tsv PR_Wagle_ASC_Pairs_Capture_All_Pairs.aggregated_case_sample.seg 2.Gistic2_Analysis_ReCapSeg - sh /xchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2010a /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Gistic2_Analysis/broadinstitute.org/cancer.genome.analysis/00264/129 /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Gistic2_Analysis/broadinstitute.org/cancer.genome.analysis/00264/129 /call_gistic2_4 /fh/subscription- XRC/tsvCat/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33645478/PR_Wagle_ASC_Pairs_Capture_All_Pairs.aggregated_case_sample.seg /xchip/gistic/variables/WES_CRSP_ICE_hg19_wex_illumina_v1_naX_Y_MT1.markers.txt /xchip/gistic/variables/hg19/hg19_with_miR_20120227.mat /xchip/gistic/CNV/blood_normals/CNV.hg19_111204/CNV.hg19.bypots.111213.txt 0.3 0.3 2 0.5 1 0.99 10 1 10000 1 mean /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/Gistic2_Analysis/broadinstitute.org/cancer.genome.analysis/00264/129 /version.txt 3.Gistic2_Report_ReCapSeg - java -Dr_flags=--vanilla -cp /xchip/tcga/gdac_prod/applications/genepattern/RunR/ RunR /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GDAC_Gistic2Report/broadinstitute.org/cancer.genome.analysis/0026 9/63/Gistic2Report.R main -L /xchip/tcga/Tools/Nozzle/v1.current - g /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GDAC_Gistic2Report/broadinstitute.org/cancer.genome.analysis/002 69/63/gene.list.txt -a /fh/subscription-XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/amp_genes.conf.99.txt - d /fh/subscription-XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/del_genes.conf.99.txt -A /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/amp_qplot.png -D /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/del_qplot.png -S /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/raw_copy_number.png -V /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/gisticVersion.txt -T /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/arraylistfile.txt -b /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/broad_significance_results.txt -l /fh/subscription- XRC/Gistic2_Analysis/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33646697/gisticInputs.txt </pre>
CoMut	Mutation_Si	<pre> 1.MutSigRun2CV- /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MutSig_2CV/broadinstitute.org/cancer.genome.analysis/10594/36/cho ose_dotkit Python-2.7 </pre>

	gnificance_2 CV	<pre> /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MutSig_2CV/broadinstitute.org/cancer.genome.analysis/10594/36/MutSig_Wrapper.sh MutSig_2CV_standalone_v3.1 /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MutSig_2CV/broadinstitute.org/cancer.genome.analysis/10594/36/Merge_MAFs.py --maf1 /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/maf1.input.tsv --maf2 /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/maf2.input.tsv PR_Wagle_ASC_Pairs_Capture_All_Pairs sh /xchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2013a /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/MutSig_2CV/broadinstitute.org/cancer.genome.analysis/10594/36/MutSig_2CV_standalone_v3.1 /xchip/cga/reference/mutsig_params/pancan_mutation_blacklist.v14 hg19.add_ERpos.CCPM.final.txt /xchip/cga/reference/mutsig_params/mutation_type_dictionary.v4.txt /xchip/cga/reference/mutsig_params/coverage_models.v5a.mat /xchip/cga/reference/mutsig_params/coverage_basewise.fwb /xchip/cga/reference/mutsig_params/target_list.hg19.v1a.txt /xchip/cga/reference/mutsig_params/context_and_effect.c65e29b.fwb /xchip/cga/reference/mutsig_params/context_and_effect.c65e29b.txt /xchip/cga/reference/mutsig_params/covariates_transformed.v5a.txt /xchip/cga/reference/mutsig_params/conservation46.fwb /xchip/cga/reference/mutsig_params/FixedWidthBinary.jar hg19 /xchip/cga/reference/mutsig_params/MutSig2CV.hg19.params.txt 2. GenerateStickFiguresMutSig2CV - /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GenerateStickFigures/broadinstitute.org/cancer.genome.analysis/00497/16/choose_dotkit Python-2.7 /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GenerateStickFigures/broadinstitute.org/cancer.genome.analysis/00497/16/generate_figs.py -s 0.99 -e True /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/PR_Wagle_ASC_Pairs_Capture_All_Pairs.final_analysis_set.maf /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/sig_genes.txt /xchip/cga/reference/annotation/db/uniprot/Uniprot_SWISS_human_Release_2010_12.shelf sh /xchip/tcga/gdac_prod/source/analysis_pipeline/scripts/choose_run_matlab2.sh -v matlab-2012a /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GenerateStickFigures/broadinstitute.org/cancer.genome.analysis/00497/16/ /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GenerateStickFigures/broadinstitute.org/cancer.genome.analysis/00497/16/mutfig 3. CoMut2CV - /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CoMut/broadinstitute.org/cancer.genome.analysis/00488/17/choose_dotkit R-2.15 Rscript /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CoMut/broadinstitute.org/cancer.genome.analysis/00488/17/coMut.R -v -o - --firehose.mode --png --sort.by.mutation.status -a PR_Wagle_ASC_Pairs_Capture_All_Pairs -s /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/sig_genes.txt -c /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/patient_counts_and_rates.txt -m /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/PR_Wagle_ASC_Pairs_Capture_All_Pairs.final_analysis_set.maf --firehose.mutsig.mutcat /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/mutcatcats.txt -q 0.75 4. MutSigNozzleReport2CV - /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GDAC_MutSig_Report/broadinstitute.org/cancer.genome.analysis/00412/47/choose_dotkit R-2.15 Rscript --slave --vanilla /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GDAC_MutSig_Report/broadinstitute.org/cancer.genome.analysis/00412/47/MutSig_Report.R -k /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/patient_counts_and_rates.txt -l /fh/subscription-XRC/CoMut/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33640421/PR_Wagle_ASC_Pairs_Capture_All_Pairs.coMut.png -m /fh/subscription-XRC/CoMut/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33640421/PR_Wagle_ASC_Pairs_Capture_All_Pairs.coMut.pdf -n /fh/subscription-XRC/GenerateStickFigures/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33640314/generatedStickFigurePaths.txt -o /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/PR_Wagle_ASC_Pairs_Capture_All_Pairs.final_analysis_set.maf /xchip/tcga/Tools/Nozzle/v1.current PR_Wagle_ASC_Pairs_Capture_All_Pairs /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/sig_genes.txt /fh/subscription-XRC/MutSig_2CV/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639571/MutSig_version.txt </pre>
Germline SNV	Germline Pipeline Simplified	<pre> 1.HaplotypeCallerSingleSampleGVCF for Pairs on Capture - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/HaplotypeCallerSingleSampleGVCF/broadinstitute.org/cancer.genome.analysis/10614/9/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/HaplotypeCallerSingleSampleGVCF/broadinstitute.org/cancer.genome.analysis/10614/9/ /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam /fh/subscription-XRC/HaplotypeCallerSingleSampleGVCF/ASC-ASCProject_Tumor_Normal/27003021/tmp ASC-ASCProject_Normal_Sample 10 /xchip/cga/reference/hg19/whole_exome_agilent.1.1_refseq_plus_3_boosters_plus_10bp_padding_minus_mito.Homo_sapiens_assembly19.targets.interval_list /xchip/cga/reference/hg19/dbsnp.134_b37.leftAligned.vcf /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta 2.Combine and Genotype GVCFs - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CombineGenotypeGVCFs/broadinstitute.org/cancer.genome.analysis/10674/1/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CombineGenotypeGVCFs/broadinstitute.org/cancer.genome.analysis/10674/1/ /fh/subscription-XRC/CombineGenotypeGVCFs/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33627528/gvcs.list.input.tsv 300 PR_Wagle_ASC_Pairs_Capture_All_Pairs 20 /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta TRUE 3.VariantRecalibrator for Germline SNPs- sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKVariantRecalibrator/broadinstitute.org/cancer.genome.analysis/10677/3/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKVariantRecalibrator/broadinstitute.org/cancer.genome.analysis/10677/3/ /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta /fh/subscription-XRC/GATKVariantRecalibrator/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33638806/tmp /fh/subscription-XRC/CombineGenotypeGVCFs/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33627528/PR_Wagle_ASC_Pairs_Capture_All_Pairs.FinalMergedGVC F.gvcf.gz /xchip/cga/reference/hg19/whole_exome_agilent.1.1_refseq_plus_3_boosters_plus_10bp_padding_minus_mito.Homo_sapiens_assembly19.targets.interval_list PR_Wagle_ASC_Pairs_Capture_All_Pairs.SNP 4.ApplyRecalibration for Germline SNPs - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKApplyRecalibration/broadinstitute.org/cancer.genome.analysis/10678/1/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKApplyRecalibration/broadinstitute.org/cancer.genome.analysis/10678/1/ /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta /fh/subscription-XRC/GATKApplyRecalibration/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639082/tmp /fh/subscription-XRC/CombineGenotypeGVCFs/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33627528/PR_Wagle_ASC_Pairs_Capture_All_Pairs.FinalMergedGVC F.gvcf.gz 99.0 /fh/subscription-XRC/GATKVariantRecalibrator/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33638806/PR_Wagle_ASC_Pairs_Capture_All_Pairs.SNP.tranches /fh/subscription-XRC/GATKVariantRecalibrator/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33638806/PR_Wagle_ASC_Pairs_Capture_All_Pairs.SNP.recal PR_Wagle_ASC_Pairs_Capture_All_Pairs.SNP 5.VariantRecalibrator for Germline Indels - </pre>

		<pre> sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKVariantRecalibrator/broadinstitute.org/cancer.genome.analysis/10677/4/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKVariantRecalibrator/broadinstitute.org/cancer.genome.analysis/10677/4/ /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta /fh/subscription-XRC/GATKVariantRecalibrator/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639233/tmp /fh/subscription-XRC/GATKApplyRecalibration/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639082/PR_Wagle_ASC_Pairs_Capture_All_Pairs.recalibrated.vcf /xchip/cga/reference/hg19/whole_exome_agilent_1.1_refseq_plus_3_boosters_plus_10bp_padding_minus_mito.Homo_sapiens_assembly19.targets.interval_list/PR_Wagle_ASC_Pairs_Capture_All_Pairs/INDEL 6. ApplyRecalibration for Germline Indels - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKApplyRecalibration/broadinstitute.org/cancer.genome.analysis/10678/1/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GATKApplyRecalibration/broadinstitute.org/cancer.genome.analysis/10678/1/ /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta /fh/subscription-XRC/GATKApplyRecalibration/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639324/tmp /fh/subscription-XRC/GATKApplyRecalibration/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639082/PR_Wagle_ASC_Pairs_Capture_All_Pairs.recalibrated.vcf 99.0 /fh/subscription-XRC/GATKVariantRecalibrator/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639233/PR_Wagle_ASC_Pairs_Capture_All_Pairs.INDEL.tranches /fh/subscription-XRC/GATKVariantRecalibrator/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639233/PR_Wagle_ASC_Pairs_Capture_All_Pairs.INDEL.recal /PR_Wagle_ASC_Pairs_Capture_All_Pairs/INDEL 7.Subset Filtered VCF to Original Set sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GermlinePipelineSubsetVCF/broadinstitute.org/cancer.genome.analysis/10683/5/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GermlinePipelineSubsetVCF/broadinstitute.org/cancer.genome.analysis/10683/5/ /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta /fh/subscription-XRC/GermlinePipelineSubsetVCF/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639325/squid_id_list.input.tsv /fh/subscription-XRC/GATKApplyRecalibration/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639324/PR_Wagle_ASC_Pairs_Capture_All_Pairs.recalibrated.vcf /PR_Wagle_ASC_Pairs_Capture_All_Pairs 8.Subset Filtered VCF to Single Sample - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GermlinePipelineSubsetVCFSingleSample/broadinstitute.org/cancer.genome.analysis/10684/10/script.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/GermlinePipelineSubsetVCFSingleSample/broadinstitute.org/cancer.genome.analysis/10684/10/ /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta /fh/subscription-XRC/GermlinePipelineSubsetVCFSingleSample/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639326/list_of_squid_ids.input.tsv /fh/subscription-XRC/GermlinePipelineSubsetVCF/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639325/PR_Wagle_ASC_Pairs_Capture_All_Pairs.vcf 9.Single Sample VCFtoMAF - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/VCFtoMAF/broadinstitute.org/cancer.genome.analysis/10649/3/vcf_to_maf.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/VCFtoMAF/broadinstitute.org/cancer.genome.analysis/10649/3/ /fh/subscription-XRC/GermlinePipelineSubsetVCFSingleSample/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639326/ASCProject_pt_id_SALIVA.GATK_PrimitiveAlleles.vcf/ASC-ASCProject_Normal_Sample/TRUE </pre>
Allelic CapSeg	AllelicCapseg	<pre> 1.VCFtoHetIntervals - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/VCFtoHetIntervals/broadinstitute.org/cancer.genome.analysis/10733/5/VCFtoHetIntervals.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/VCFtoHetIntervals/broadinstitute.org/cancer.genome.analysis/10733/5/ /fh/subscription-XRC/GermlinePipelineSubsetVCFSingleSample/PR_Wagle_ASC_Pairs_Capture_All_Pairs/33639326/ASCProject_pt_id_SALIVA.GATK_PrimitiveAlleles.vcf 2.HetSitePullDown - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallSomaticMutations/broadinstitute.org/cancer.genome.analysis/00004/131/runBroadJava7.sh java -Xmx2g -jar /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/CallSomaticMutations/broadinstitute.org/cancer.genome.analysis/00004/131/muTect-1.1.6.jar --analysis_type MuTect --intervals /fh/subscription-XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/33640808/iteration1/ /gatk-scatter.0000000015.interval_list --normal_sample_name ASC-ASCProject_Normal_Sample -I-normal /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_SALIVA/v2/ASCProject_pt_id_SALIVA.bam --tumor_sample_name ASC-ASCProject_Tumor_Sample-I-tumor /seq/picard_aggregation/RP-1156/Exome/ASCProject_pt_id_T1/v2/ASCProject_pt_id_T1.bam --reference_sequence /seq/references/Homo_sapiens_assembly19/v1/Homo_sapiens_assembly19.fasta --dbsnp /xchip/cga/reference/hg19/dbsnp_134_b37.leftAligned.vcf --cosmic /xchip/cga/reference/hg19/hg19_cosmic_v54_120711.vcf --normal_panel /xchip/cga/reference/hg19/refseq_exome_10bp_hg19_300_1kg_normal_panel.vcf --out ASC-ASCProject_Tumor_Normal.call_stats.txt --coverage_file ASC-ASCProject_Tumor_Normal.coverage.wig.txt --power_file ASC-ASCProject_Tumor_Normal.power.wig.txt --downsample_to_coverage 100000 --enable_extended_output --fraction_contamination 0.007 --force_output 3. HetPullDownPostProcess - sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/HetPullDownPostprocess/broadinstitute.org/cancer.genome.analysis/10621/10/CallStatsToCov.sh /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/HetPullDownPostprocess/broadinstitute.org/cancer.genome.analysis/10621/10/ /fh/subscription-XRC/CallSomaticMutations/ASC-ASCProject_Tumor_Normal/33640808/iteration1/ASC-ASCProject_Tumor_Normal.call_stats.txt 4. AllelicCapseg - /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/AllelicCapseg/broadinstitute.org/cancer.genome.analysis/10622/38/choose_dotkit R-2.15 Rscript /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/AllelicCapseg/broadinstitute.org/cancer.genome.analysis/10622/38/AllelicCapseg_cli.R --SID=ASC-ASCProject_Tumor_Sample --capseg.probe.fn=/fh/subscription-XRC/recapseg_tumor_pcov/ASC-ASCProject_Tumor_Sample/27011228/recapseg.tn.ASC-ASCProject_Tumor_Sample.tsv --capseg.seg.fn=/fh/subscription-XRC/recapseg_tumor_pcov/ASC-ASCProject_Tumor_Sample/27011228/recapseg.ASC-ASCProject_Tumor_Sample.seg --germline.het.fn=/fh/subscription-XRC/HetPullDownPostprocess/ASC-ASCProject_Tumor_Normal/33641000/Tumor.cov --drop.x=FALSE --drop.y=TRUE --seg.merge.thresh=0.5 --min.seg.size=3 --verbose=TRUE --base.output.dir= --initial.merge=TRUE --split.merge=TRUE --outlier.thresh=0.005 --working.dir=/xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/AllelicCapseg/broadinstitute.org/cancer.genome.analysis/10622/38/ </pre>

ABSOLUTE	ABSOLUTE _v1.5_WES	<pre>/broad/software/free/Linux/redhat_5_x86_64/pkgs/r_2.15.3/bin/Rscript /xchip/tcga/gdac_prod/applications/process_mgmt/firehose_task_registry/cga/ABSOLUTE_v1.5/broadinstitute.org/cancer.genome.analysis/10971/2 3//ABSOLUTE_cli_start.R --seg_dat_fn /fh/subscription-XRC/AllelicCapseg/ASC-ASCProject_Tumor_Normal/33641022/results/ASC- ASCProject_Tumor_Sample.tsv --maf_fn /fh/subscription-XRC/PoN_filter_py/ASC-ASCProject_Tumor_Normal/30688317/ASC- ASCProject_Tumor_Normal.noalign.PoNfiltered_agilent_ICE.maf --indelmaf_fn /fh/subscription-XRC/MafMasterFilter/ASC- ASCProject_Tumor_Normal/31668062/ASC-ASCProject_Tumor_Normal.pon_filtered.txt --sample_name ASC-ASCProject_Tumor_Normal -- results_dir raw_results/ --ssnv_skew 0.9568517 --abs_lib_dir /xchip/tcga/Tools/absolute/releases/v1.5/</pre>
----------	-------------------------------	--