





Saliva collection instructions


Do NOT eat, drink, smoke or chew gum for 30 minutes before giving your saliva sample.
Do NOT remove the plastic film from the funnel lid.




1. Spit until the amount of saliva* (not bubbles) reaches the fill line.




2. Close lid tightly by pushing down hard on the funnel lid until you hear a loud click.



3. Hold the tube upright. Unscrew the funnel from the tube.



4. Use the small cap to close the tube tightly.



5. Shake the capped tube for 5 seconds. Discard funnel.

⚠ Small cap, choking hazard. Wash with water if stabilizing liquid comes in contact with eyes or skin. Do NOT ingest.

Mailing instructions



6. Locate specimen bag provided. Do NOT remove absorbent pad.



7. Seal the capped tube into the bag.



8. Place bag with sample back into original box.



9. Close box and seal shut.



10. Mail from nearest postal location.



Need help?

Phone: 651-293-5029

Email: info@mpcproject.org



Metastatic
Prostate Cancer
Project

www.mpcproject.org

PD-PR-00648 Issue 2/2017-09